



HydroMF Fitting Guide

Fitting Procedure

- Full refraction and eye examination.
- Keratometry.
- Select diameter 2.00-3.00mm larger than HVID. 14.00mm or 14.50mm available.
- Select the closest base curve available based on K readings:

| K readings (mm) | Base Curve |
|-----------------|------------|
| 7.45 and below | 8.30 |
| 7.45 to 7.90 | 8.70 |
| 7.90 and above | 9.10 |

- Determine the dominant eye.
- Select power from PV conversion tables dominant (Distance) non-dominant (near); Distance power = Best Vision Sphere with BVD adjustment.
- Order as PV 1.50D, PV 2.00D or PV 2.50D Design (as per table).

Lens Fit

- Insert lens and allow to settle for 10 minutes.

| Characteristic | Good Fit | Steep Fit | Flat Fit |
|---------------------|--|--|---------------------------------|
| Comfort | Good | Good initially | Poor |
| Centration | Good | Often good | Often poor |
| Coverage | Full limbal coverage in all directions of gaze | Often good | Poor |
| Edge | Good | May cause conjunctival indentation/blanching | Edge stand-off/fluting |
| Movement with blink | Good | Inadequate- None | Excessive |
| Push up test | Easily moved/recentres | Difficult to displace/slow recenteration | Easily moved/poor recenteration |

Continue to the following page to see the next steps.

Good fit - Primary Position



Good fit - Nasal Position



Steep fit - Primary Position



Steep fit - Nasal Position



Flat fit - Primary Position



Flat fit - Nasal Position



- Check binocular acuity for near first, then intermediate and distance in normal light without occluding either eye.
- If distance and near binocular acuities are within a line or two of the patient's potential, leave the lenses in for another 45 to 50 minutes to equilibrate the eyes.
- If correction is needed, binocularly over-refract the reading (near) eye first, then the distance eye without occluding either eye or using a phoropter.
- Order an exchange lens if required.

Tips

- Use a best vision sphere - transpose cylinders into minus form then add half of the cylinder to the spherical distance correction before selecting power from the chart.
- Because of the 20.00D power spread, 0.25D can have unusually positive effects on acuity.
- When it is difficult to reach a satisfactory balance between distance and vision, try a flatter lens.
- If the patient is emmetropic (no distance correction is required) the patient may be satisfied with one lens in their non-dominant eye.
- The lens will be labelled with the combined power form (i.e. Distance Power + Reading Addition).

Modality

HydroMF lenses are designed to last up to 12 months. The lenses may need to be replaced more frequently due to the patient's wear and care; this is at the discretion of the practitioner.

Care

Hydrogen peroxide systems and multipurpose systems may be used. Please contact us for information on caring for lenses.

All measurements are in Dioptres.

PV Multifocal Lenses – Dominant Eye (Distance)

| Spectacle Add | +1.25 | +1.50 | +1.75 | +2.00 | +2.25 | +2.50 | +2.75 | |
|---------------|--------------|-------|--------------|-------|-------|--------------|-------|-------|
| Dist | +1.50 PV ADD | | +2.00 PV ADD | | | +2.50 PV ADD | | |
| 0.00 | +1.00 | +1.25 | +1.50 | +1.75 | +1.75 | +2.25 | +2.25 | +2.50 |
| -0.25 | +0.75 | +1.00 | +1.25 | +1.50 | +1.50 | +2.00 | +2.00 | +2.25 |
| -0.50 | +0.50 | +0.75 | +1.00 | +1.25 | +1.25 | +1.75 | +1.75 | +2.00 |
| -0.75 | +0.25 | +0.50 | +0.75 | +1.00 | +1.00 | +1.50 | +1.50 | +1.75 |
| -1.00 | 0.00 | +0.25 | +0.50 | +0.75 | +0.75 | +1.25 | +1.25 | +1.50 |
| -1.25 | -0.25 | 0.00 | +0.25 | +0.50 | +0.50 | +1.00 | +1.00 | +1.25 |
| -1.50 | -0.50 | -0.25 | 0.00 | +0.25 | +0.25 | +0.75 | +0.75 | +1.00 |
| -1.75 | -0.75 | -0.50 | -0.25 | 0.00 | 0.00 | +0.50 | +0.50 | +0.75 |
| -2.00 | -1.00 | -0.75 | -0.50 | -0.25 | -0.25 | +0.25 | +0.25 | +0.50 |
| -2.25 | -1.25 | -1.00 | -0.75 | -0.50 | -0.50 | 0.00 | 0.00 | +0.25 |
| -2.50 | -1.50 | -1.25 | -1.00 | -0.75 | -0.75 | -0.25 | -0.25 | 0.00 |
| -2.75 | -1.75 | -1.50 | -1.25 | -1.00 | -1.00 | -0.50 | -0.50 | -0.25 |
| -3.00 | -2.00 | -1.75 | -1.50 | -1.25 | -1.25 | -0.75 | -0.75 | -0.50 |
| -3.25 | -2.25 | -2.00 | -1.75 | -1.50 | -1.50 | -1.00 | -1.00 | -0.75 |
| -3.50 | -2.50 | -2.25 | -2.00 | -1.75 | -1.75 | -1.25 | -1.25 | -1.00 |
| -3.75 | -2.75 | -2.50 | -2.25 | -2.00 | -2.00 | -1.50 | -1.50 | -1.25 |
| -4.00 | -3.00 | -2.75 | -2.50 | -2.25 | -2.25 | -1.75 | -1.75 | -1.50 |
| -4.25 | -3.25 | -3.00 | -2.75 | -2.50 | -2.50 | -2.00 | -2.00 | -1.75 |
| -4.50 | -3.50 | -3.25 | -3.00 | -2.75 | -2.75 | -2.25 | -2.25 | -2.00 |
| -4.75 | -3.75 | -3.50 | -3.25 | -3.00 | -3.00 | -2.50 | -2.50 | -2.25 |
| -5.00 | -4.00 | -3.75 | -3.50 | -3.25 | -3.25 | -2.75 | -2.75 | -2.50 |
| -5.25 | -4.25 | -4.00 | -3.75 | -3.50 | -3.50 | -3.00 | -3.00 | -2.75 |
| -5.50 | -4.50 | -4.25 | -4.00 | -3.75 | -3.75 | -3.25 | -3.25 | -3.00 |
| -5.75 | -4.75 | -4.50 | -4.25 | -4.00 | -4.00 | -3.50 | -3.50 | -3.25 |
| -6.00 | -5.00 | -4.75 | -4.50 | -4.25 | -4.25 | -3.75 | -3.75 | -3.50 |
| -6.25 | -5.25 | -5.00 | -4.75 | -4.50 | -4.50 | -4.00 | -4.00 | -3.75 |
| -6.50 | -5.50 | -5.25 | -5.00 | -4.75 | -4.75 | -4.25 | -4.25 | -4.00 |
| -6.75 | -5.75 | -5.50 | -5.25 | -5.00 | -5.00 | -4.50 | -4.50 | -4.25 |
| -7.00 | -6.00 | -5.75 | -5.50 | -5.25 | -5.25 | -4.75 | -4.75 | -4.50 |
| -7.25 | -6.25 | -6.00 | -5.75 | -5.50 | -5.50 | -5.00 | -5.00 | -4.75 |
| -7.50 | -6.50 | -6.25 | -6.00 | -5.75 | -5.75 | -5.25 | -5.25 | -5.00 |
| -7.75 | -6.75 | -6.50 | -6.25 | -6.00 | -6.00 | -5.50 | -5.50 | -5.25 |
| -8.00 | -7.00 | -6.75 | -6.50 | -6.25 | -6.25 | -5.75 | -5.75 | -5.50 |
| -8.25 | -7.25 | -7.00 | -6.75 | -6.50 | -6.50 | -6.00 | -6.00 | -5.75 |
| -8.50 | -7.50 | -7.25 | -7.00 | -6.75 | -6.75 | -6.25 | -6.25 | -6.00 |
| -8.75 | -7.75 | -7.50 | -7.25 | -7.00 | -7.00 | -6.50 | -6.50 | -6.25 |
| -9.00 | -8.00 | -7.75 | -7.50 | -7.25 | -7.25 | -6.75 | -6.75 | -6.50 |
| -9.25 | -8.25 | -8.00 | -7.75 | -7.50 | -7.50 | -7.00 | -7.00 | -6.75 |
| -9.50 | -8.50 | -8.25 | -8.00 | -7.75 | -7.75 | -7.25 | -7.25 | -7.00 |
| -9.75 | -8.75 | -8.50 | -8.25 | -8.00 | -8.00 | -7.50 | -7.50 | -7.25 |
| -10.00 | -9.00 | -8.75 | -8.50 | -8.25 | -8.25 | -7.75 | -7.75 | -7.50 |

PV Multifocal Lenses – Non-Dominant Eye (Near)

| Add | +1.25 | +1.50 | +1.75 | +2.00 | +2.25 | +2.50 | +2.75 | |
|--------|--------------|-------|--------------|-------|-------|--------------|-------|-------|
| Dist | +1.50 PV ADD | | +2.00 PV ADD | | | +2.50 PV ADD | | |
| 0.00 | +1.25 | +1.50 | +1.75 | +2.00 | +2.00 | +2.50 | +2.50 | +2.75 |
| -0.25 | +1.00 | +1.25 | +1.50 | +1.75 | +1.75 | +2.25 | +2.25 | +2.50 |
| -0.50 | +0.75 | +1.00 | +1.25 | +1.50 | +1.50 | +2.00 | +2.00 | +2.25 |
| -0.75 | +0.50 | +0.75 | +1.00 | +1.25 | +1.25 | +1.75 | +1.75 | +2.00 |
| -1.00 | +0.25 | +0.50 | +0.75 | +1.00 | +1.00 | +1.50 | +1.50 | +1.75 |
| -1.25 | 0.00 | +0.25 | +0.50 | +0.75 | +0.75 | +1.25 | +1.25 | +1.50 |
| -1.50 | -0.25 | 0.00 | +0.25 | +0.50 | +0.50 | +1.00 | +1.00 | +1.25 |
| -1.75 | -0.50 | -0.25 | 0.00 | +0.25 | +0.25 | +0.75 | +0.75 | +1.00 |
| -2.00 | -0.75 | -0.50 | -0.25 | 0.00 | 0.00 | +0.50 | +0.50 | +0.75 |
| -2.25 | -1.00 | -0.75 | -0.50 | -0.25 | -0.25 | +0.25 | +0.25 | +0.50 |
| -2.50 | -1.25 | -1.00 | -0.75 | -0.50 | -0.50 | 0.00 | 0.00 | +0.25 |
| -2.75 | -1.50 | -1.25 | -1.00 | -0.75 | -0.75 | -0.25 | -0.25 | 0.00 |
| -3.00 | -1.75 | -1.50 | -1.25 | -1.00 | -1.00 | -0.50 | -0.50 | -0.25 |
| -3.25 | -2.00 | -1.75 | -1.50 | -1.25 | -1.25 | -0.75 | -0.75 | -0.50 |
| -3.50 | -2.25 | -2.00 | -1.75 | -1.50 | -1.50 | -1.00 | -1.00 | -0.75 |
| -3.75 | -2.50 | -2.25 | -2.00 | -1.75 | -1.75 | -1.25 | -1.25 | -1.00 |
| -4.00 | -2.75 | -2.50 | -2.25 | -2.00 | -2.00 | -1.50 | -1.50 | -1.25 |
| -4.25 | -3.00 | -2.75 | -2.50 | -2.25 | -2.25 | -1.75 | -1.75 | -1.50 |
| -4.50 | -3.25 | -3.00 | -2.75 | -2.50 | -2.50 | -2.00 | -2.00 | -1.75 |
| -4.75 | -3.50 | -3.25 | -3.00 | -2.75 | -2.75 | -2.25 | -2.25 | -2.00 |
| -5.00 | -3.75 | -3.50 | -3.25 | -3.00 | -3.00 | -2.50 | -2.50 | -2.25 |
| -5.25 | -4.00 | -3.75 | -3.50 | -3.25 | -3.25 | -2.75 | -2.75 | -2.50 |
| -5.50 | -4.25 | -4.00 | -3.75 | -3.50 | -3.50 | -3.00 | -3.00 | -2.75 |
| -5.75 | -4.50 | -4.25 | -4.00 | -3.75 | -3.75 | -3.25 | -3.25 | -3.00 |
| -6.00 | -4.75 | -4.50 | -4.25 | -4.00 | -4.00 | -3.50 | -3.50 | -3.25 |
| -6.25 | -5.00 | -4.75 | -4.50 | -4.25 | -4.25 | -3.75 | -3.75 | -3.50 |
| -6.50 | -5.25 | -5.00 | -4.75 | -4.50 | -4.50 | -4.00 | -4.00 | -3.75 |
| -6.75 | -5.50 | -5.25 | -5.00 | -4.75 | -4.75 | -4.25 | -4.25 | -4.00 |
| -7.00 | -5.75 | -5.50 | -5.25 | -5.00 | -5.00 | -4.50 | -4.50 | -4.25 |
| -7.25 | -6.00 | -5.75 | -5.50 | -5.25 | -5.25 | -4.75 | -4.75 | -4.50 |
| -7.50 | -6.25 | -6.00 | -5.75 | -5.50 | -5.50 | -5.00 | -5.00 | -4.75 |
| -7.75 | -6.50 | -6.25 | -6.00 | -5.75 | -5.75 | -5.25 | -5.25 | -5.00 |
| -8.00 | -6.75 | -6.50 | -6.25 | -6.00 | -6.00 | -5.50 | -5.50 | -5.25 |
| -8.25 | -7.00 | -6.75 | -6.50 | -6.25 | -6.25 | -5.75 | -5.75 | -5.50 |
| -8.50 | -7.25 | -7.00 | -6.75 | -6.50 | -6.50 | -6.00 | -6.00 | -5.75 |
| -8.75 | -7.50 | -7.25 | -7.00 | -6.75 | -6.75 | -6.25 | -6.25 | -6.00 |
| -9.00 | -7.75 | -7.50 | -7.25 | -7.00 | -7.00 | -6.50 | -6.50 | -6.25 |
| -9.25 | -8.00 | -7.75 | -7.50 | -7.25 | -7.25 | -6.75 | -6.75 | -6.50 |
| -9.50 | -8.25 | -8.00 | -7.75 | -7.50 | -7.50 | -7.00 | -7.00 | -6.75 |
| -9.75 | -8.50 | -8.25 | -8.00 | -7.75 | -7.75 | -7.25 | -7.25 | -7.00 |
| -10.00 | -8.75 | -8.50 | -8.25 | -8.00 | -8.00 | -7.50 | -7.50 | -7.25 |

PV Multifocal Lenses –Dominant Eye (Distance)

| Add | +1.25 | +1.50 | +1.75 | +2.00 | +2.25 | +2.50 | +2.75 | |
|--------|--------------|--------|--------------|--------|--------|--------------|--------|--------|
| Dist | +1.50 PV ADD | | +2.00 PV ADD | | | +2.50 PV ADD | | |
| +0.25 | +1.50 | +1.50 | +2.00 | +2.00 | +2.00 | +2.50 | +2.50 | +2.75 |
| +0.50 | +1.75 | +1.75 | +2.25 | +2.25 | +2.25 | +2.75 | +2.75 | +3.00 |
| +0.75 | +2.00 | +2.00 | +2.50 | +2.50 | +2.50 | +3.00 | +3.00 | +3.25 |
| +1.00 | +2.25 | +2.25 | +2.75 | +2.75 | +2.75 | +3.25 | +3.25 | +3.50 |
| +1.25 | +2.50 | +2.50 | +3.00 | +3.00 | +3.00 | +3.50 | +3.50 | +3.75 |
| +1.50 | +2.75 | +2.75 | +3.25 | +3.25 | +3.25 | +3.75 | +3.75 | +4.00 |
| +1.75 | +3.00 | +3.00 | +3.50 | +3.50 | +3.50 | +4.00 | +4.00 | +4.25 |
| +2.00 | +3.25 | +3.25 | +3.75 | +3.75 | +3.75 | +4.25 | +4.25 | +4.50 |
| +2.25 | +3.50 | +3.50 | +4.00 | +4.00 | +4.00 | +4.50 | +4.50 | +4.75 |
| +2.50 | +3.75 | +3.75 | +4.25 | +4.25 | +4.25 | +4.75 | +4.75 | +5.00 |
| +2.75 | +4.00 | +4.00 | +4.50 | +4.50 | +4.50 | +5.00 | +5.00 | +5.25 |
| +3.00 | +4.25 | +4.25 | +4.75 | +4.75 | +4.75 | +5.25 | +5.25 | +5.50 |
| +3.25 | +4.50 | +4.50 | +5.00 | +5.00 | +5.00 | +5.50 | +5.50 | +5.75 |
| +3.50 | +4.75 | +4.75 | +5.25 | +5.25 | +5.25 | +5.75 | +5.75 | +6.00 |
| +3.75 | +5.00 | +5.00 | +5.50 | +5.50 | +5.50 | +6.00 | +6.00 | +6.25 |
| +4.00 | +5.25 | +5.25 | +5.75 | +5.75 | +5.75 | +6.25 | +6.25 | +6.50 |
| +4.25 | +5.50 | +5.50 | +6.00 | +6.00 | +6.00 | +6.50 | +6.50 | +6.75 |
| +4.50 | +5.75 | +5.75 | +6.25 | +6.25 | +6.25 | +6.75 | +6.75 | +7.00 |
| +4.75 | +6.00 | +6.00 | +6.50 | +6.50 | +6.50 | +7.00 | +7.00 | +7.25 |
| +5.00 | +6.25 | +6.25 | +6.75 | +6.75 | +6.75 | +7.25 | +7.25 | +7.50 |
| +5.25 | +6.50 | +6.50 | +7.00 | +7.00 | +7.00 | +7.50 | +7.50 | +7.75 |
| +5.50 | +6.75 | +6.75 | +7.25 | +7.25 | +7.25 | +7.75 | +7.75 | +8.00 |
| +5.75 | +7.00 | +7.00 | +7.50 | +7.50 | +7.50 | +8.00 | +8.00 | +8.25 |
| +6.00 | +7.25 | +7.25 | +7.75 | +7.75 | +7.75 | +8.25 | +8.25 | +8.50 |
| +6.25 | +7.50 | +7.50 | +8.00 | +8.00 | +8.00 | +8.50 | +8.50 | +8.75 |
| +6.50 | +7.75 | +7.75 | +8.25 | +8.25 | +8.25 | +8.75 | +8.75 | +9.00 |
| +6.75 | +8.00 | +8.00 | +8.50 | +8.50 | +8.50 | +9.00 | +9.00 | +9.25 |
| +7.00 | +8.25 | +8.25 | +8.75 | +8.75 | +8.75 | +9.25 | +9.25 | +9.50 |
| +7.25 | +8.50 | +8.50 | +9.00 | +9.00 | +9.00 | +9.50 | +9.50 | +9.75 |
| +7.50 | +8.75 | +8.75 | +9.25 | +9.25 | +9.25 | +9.75 | +9.75 | +10.00 |
| +7.75 | +9.00 | +9.00 | +9.50 | +9.50 | +9.50 | +10.00 | +10.00 | +10.25 |
| +8.00 | +9.25 | +9.25 | +9.75 | +9.75 | +9.75 | +10.25 | +10.25 | +10.50 |
| +8.25 | +9.50 | +9.50 | +10.00 | +10.00 | +10.00 | +10.50 | +10.50 | +10.75 |
| +8.50 | +9.75 | +9.75 | +10.25 | +10.25 | +10.25 | +10.75 | +10.75 | +11.00 |
| +8.75 | +10.00 | +10.00 | +10.50 | +10.50 | +10.50 | +11.00 | +11.00 | +11.25 |
| +9.00 | +10.25 | +10.25 | +10.75 | +10.75 | +10.75 | +11.25 | +11.25 | +11.50 |
| +9.25 | +10.50 | +10.50 | +11.00 | +11.00 | +11.00 | +11.50 | +11.50 | +11.75 |
| +9.50 | +10.75 | +10.75 | +11.25 | +11.25 | +11.25 | +11.75 | +11.75 | +12.00 |
| +9.75 | +11.00 | +11.00 | +11.50 | +11.50 | +11.50 | +12.00 | +12.00 | +12.25 |
| +10.00 | +11.25 | +11.25 | +11.75 | +11.75 | +11.75 | +12.25 | +12.25 | +12.50 |

PV Multifocal Lenses – Non-Dominant Eye (Near)

| Add | +1.25 | +1.50 | +1.75 | +2.00 | +2.25 | +2.50 | +2.75 | |
|--------|--------------|--------|--------------|--------|--------|--------------|--------|--------|
| Dist | +1.50 PV ADD | | +2.00 PV ADD | | | +2.50 PV ADD | | |
| +0.25 | +1.75 | +1.75 | +2.25 | +2.25 | +2.50 | +2.75 | +2.75 | +3.00 |
| +0.50 | +2.00 | +2.00 | +2.50 | +2.50 | +2.75 | +3.00 | +3.00 | +3.25 |
| +0.75 | +2.25 | +2.25 | +2.75 | +2.75 | +3.00 | +3.25 | +3.25 | +3.50 |
| +1.00 | +2.50 | +2.50 | +3.00 | +3.00 | +3.25 | +3.50 | +3.50 | +3.75 |
| +1.25 | +2.75 | +2.75 | +3.25 | +3.25 | +3.50 | +3.75 | +3.75 | +4.00 |
| +1.50 | +3.00 | +3.00 | +3.50 | +3.50 | +3.75 | +4.00 | +4.00 | +4.25 |
| +1.75 | +3.25 | +3.25 | +3.75 | +3.75 | +4.00 | +4.25 | +4.25 | +4.50 |
| +2.00 | +3.50 | +3.50 | +4.00 | +4.00 | +4.25 | +4.50 | +4.50 | +4.75 |
| +2.25 | +3.75 | +3.75 | +4.25 | +4.25 | +4.50 | +4.75 | +4.75 | +5.00 |
| +2.50 | +4.00 | +4.00 | +4.50 | +4.50 | +4.75 | +5.00 | +5.00 | +5.25 |
| +2.75 | +4.25 | +4.25 | +4.75 | +4.75 | +5.00 | +5.25 | +5.25 | +5.50 |
| +3.00 | +4.50 | +4.50 | +5.00 | +5.00 | +5.25 | +5.50 | +5.50 | +5.75 |
| +3.25 | +4.75 | +4.75 | +5.25 | +5.25 | +5.50 | +5.75 | +5.75 | +6.00 |
| +3.50 | +5.00 | +5.00 | +5.50 | +5.50 | +5.75 | +6.00 | +6.00 | +6.25 |
| +3.75 | +5.25 | +5.25 | +5.75 | +5.75 | +6.00 | +6.25 | +6.25 | +6.50 |
| +4.00 | +5.50 | +5.50 | +6.00 | +6.00 | +6.25 | +6.50 | +6.50 | +6.75 |
| +4.25 | +5.75 | +5.75 | +6.25 | +6.25 | +6.50 | +6.75 | +6.75 | +7.00 |
| +4.50 | +6.00 | +6.00 | +6.50 | +6.50 | +6.75 | +7.00 | +7.00 | +7.25 |
| +4.75 | +6.25 | +6.25 | +6.75 | +6.75 | +7.00 | +7.25 | +7.25 | +7.50 |
| +5.00 | +6.50 | +6.50 | +7.00 | +7.00 | +7.25 | +7.50 | +7.50 | +7.75 |
| +5.25 | +6.75 | +6.75 | +7.25 | +7.25 | +7.50 | +7.75 | +7.75 | +8.00 |
| +5.50 | +7.00 | +7.00 | +7.50 | +7.50 | +7.75 | +8.00 | +8.00 | +8.25 |
| +5.75 | +7.25 | +7.25 | +7.75 | +7.75 | +8.00 | +8.25 | +8.25 | +8.50 |
| +6.00 | +7.50 | +7.50 | +8.00 | +8.00 | +8.25 | +8.50 | +8.50 | +8.75 |
| +6.25 | +7.75 | +7.75 | +8.25 | +8.25 | +8.50 | +8.75 | +8.75 | +9.00 |
| +6.50 | +8.00 | +8.00 | +8.50 | +8.50 | +8.75 | +9.00 | +9.00 | +9.25 |
| +6.75 | +8.25 | +8.25 | +8.75 | +8.75 | +9.00 | +9.25 | +9.25 | +9.50 |
| +7.00 | +8.50 | +8.50 | +9.00 | +9.00 | +9.25 | +9.50 | +9.50 | +9.75 |
| +7.25 | +8.75 | +8.75 | +9.25 | +9.25 | +9.50 | +9.75 | +9.75 | +10.00 |
| +7.50 | +9.00 | +9.00 | +9.50 | +9.50 | +9.75 | +10.00 | +10.00 | +10.25 |
| +7.75 | +9.25 | +9.25 | +9.75 | +9.75 | +10.00 | +10.25 | +10.25 | +10.50 |
| +8.00 | +9.50 | +9.50 | +10.00 | +10.00 | +10.25 | +10.50 | +10.50 | +10.75 |
| +8.25 | +9.75 | +9.75 | +10.25 | +10.25 | +10.50 | +10.75 | +10.75 | +11.00 |
| +8.50 | +10.00 | +10.00 | +10.50 | +10.50 | +10.75 | +11.00 | +11.00 | +11.25 |
| +8.75 | +10.25 | +10.25 | +10.75 | +10.75 | +11.00 | +11.25 | +11.25 | +11.50 |
| +9.00 | +10.50 | +10.50 | +11.00 | +11.00 | +11.25 | +11.50 | +11.50 | +11.75 |
| +9.25 | +10.75 | +10.75 | +11.25 | +11.25 | +11.50 | +11.75 | +11.75 | +12.00 |
| +9.50 | +11.00 | +11.00 | +11.50 | +11.50 | +11.75 | +12.00 | +12.00 | +12.25 |
| +9.75 | +11.25 | +11.25 | +11.75 | +11.75 | +12.00 | +12.25 | +12.25 | +12.50 |
| +10.00 | +11.50 | +11.50 | +12.00 | +12.00 | +12.25 | +12.50 | +12.50 | +12.75 |