



Delta Scleral

The latest innovation in scleral design. This unique lens provides a simple diagnostic fitting process with easy-to-follow steps.

CAN AID THE FOLLOWING

- Correction of Myopia
- Correction of Hyperopia
- Correction of Astigmatism
- Correction of refractive error caused by corneal irregularities including Keratoconus

EXTRAS AVAILABLE

- Front Surface Toric
- Edge Vaults (maximum of 2)*
- Quad Haptic*

*Quad Haptics AND Edge Vaults are available to have on the same lens.

MATERIAL OPTIONS

- Optimum Extra Roflufocon D 3 (100) Blue, Clear
- Optimum Infinite Tisilfocon A 3 (180) Blue, Clear
- Fitting sets = Optimum Extra Non-UV Clear

PARAMETERS

- Diameter 15.50, 16.50, 17.50
- Power +/- 30.00D (0.25D steps)
- Cyl -15.00DC to -0.25DC (0.25DC steps)
- Axis 1° to 180°
- Central clearance, Steep Sag, Limbal Clearance, Edge alignment can all be adjusted to a maximum of +/- 1000 microns (50 micron steps).

All measurements are in mm unless otherwise stated.

MODALITY

12 months maximum. Replacement frequency is at the discretion of the eye care professional.

FITTING

A fitting set is required.

Refer to our Delta fitting guide available on our website.

DELTA 18 LENS SET

Optimum Extra Non-UV Clear

15.5			
Dx #	BC	Pwr	Sag
1	8.50	Plano	3800
2	8.10	-1.00	4100
3	7.80	-2.00	4400
4	7.50	-3.00	4700
5	7.20	-4.00	5000
6	6.95	-5.00	5300

16.5			
Dx #	BC	Pwr	Sag
7	8.70	Plano	4100
8	8.30	-1.00	4400
9	8.00	-2.00	4700
10	7.65	-3.00	5000
11	7.40	-4.00	5300
12	7.10	-5.00	5600

17.5			
Dx #	BC	Pwr	Sag
13	9.00	Plano	4400
14	8.55	-1.00	4700
15	8.15	-2.00	5000
16	7.80	-3.00	5300
17	7.50	-4.00	5600
18	7.25	-5.00	5900

