

Nissel KII RGP

A highly successful aberration controlled RGP lens for the management of Keratoconus.

CAN AID THE FOLLOWING

- Myopia
- Hyperopia
- Astigmatism
- Keratoconus

EXTRAS AVAILABLE

- Toric Flange
- Back surface toric
- Front surface toric
- Bi-Toric

DESIGN

- Spherical blended front surface
- Aspheric back surface

MATERIAL OPTIONS

- Optimum Comfort Rofluocon C 3 (65)

Limited materials available on request.

PARAMETERS

- Base Curve Steepest available 4.80 (0.05 steps)
- Diameter 8.10, 8.40, 8.70 (standard), 9.00, 9.30
- Power $-/+$ 30.00D (0.25D steps) - higher powers may be available.
- Edge lift Standard, Flat 1, Flat 2, Steep 1, Steep 2 (1 step = 0.03)
- Handling tints Blue, Clear, Green, Grey

All measurements are in mm unless otherwise stated.

FITTING

Refer to our Nissel KII RGP fitting guide available on our website.

MODALITY

12 months maximum. Replacement frequency is at the discretion of the eye care professional.

