



Zero 6 Range

A family of lathe-cut spherical and toric custom-made thin design soft lenses, suitable for everyday wear.

CAN AID THE FOLLOWING

- Myopia
- Hyperopia
- Astigmatism

DESIGN

Spherical

- Spherical blended front surface
- Spherical blended back surface

Toric

- Aspheric front surface toric with dynamic stabilisation
- Spherical blended back surface

MATERIAL OPTIONS

- Contaflex 38% water content Filcon 1 (8)
- Definitive 65% water content SiH Safrofilcon 5B (62)

FITTING

Refer to our Zero 6 fitting guide available on our website.

MODALITY

Recommended 6 months (12 months maximum).
Replacement frequency is at the discretion of the eye care professional.

PARAMETERS

Zero 6

Base Curve	8.40 to 9.30 (0.30 steps)
Diameter	14.00
Power Range Sphere	-10.00 to +10.00 D (0.25D steps) -10.00 to -20.00 D (0.50D steps) +10.00 to +20.00 D (0.50D steps)


Zero 6 Rx

Base Curve	8.10 to 9.60 (0.30 steps)
Diameter	13.50 to 15.00 (0.50 steps)
Power Range Sphere	-20.00 to +20.00D (0.25D steps) -20.00 to -30.00D (0.50D steps) +20.00 to +30.00D (0.50D steps)

Zero 6 Toric

Base Curve	8.10 to 9.60 (0.30 steps)
Diameter	14.00
Power Range Sphere	-25.00 to +25.00D (0.25D steps)
Power Range Cyl	-0.50 to -8.00D (0.25D steps)
Axis	1 to 180°

All measurements are in mm unless otherwise stated.



Toric markings on the inferior part of the lens 15° apart.
Central line should sit at 6 o'clock.

